




























# DEKEMA. Technical data.

	AUSTROMAT™ 220	AUSTROMAT™ 624	AUSTROMAT™ 624i™	AUSTROMAT™ 654	AUSTROMAT™ 654i™
Size (WxHxD)	10,47" x 20,91" x 10,16"	15,08" x 25,71" x 13,50"	15,08" x 25,71" x 13,50"	15,08" x 25,71" x 13,50"	15,08" x 25,71" x 13,50"
Weight	28,5 lbs	38,5 lbs	42 lbs	48,5 lbs	53 lbs
Power consumption	95 – 105 V/50–60 Hz or 110 – 125 V/50–60 Hz or 220 – 240 V/50–60 Hz	95 – 105 V/50–60 Hz or 110 – 125 V/50–60 Hz or 220 – 240 V/50–60 Hz	95 – 105 V/50–60 Hz or 110 – 125 V/50–60 Hz or 220 – 240 V/50–60 Hz	95 – 105 V/50–60 Hz or 110 – 125 V/50–60 Hz or 220 – 240 V/50–60 Hz	95 – 105 V/50–60 Hz or 110 – 125 V/50–60 Hz or 220 – 240 V/50–60 Hz
Maximum power consumption	1.500 W	1.500 W	1.500 W	1.500 W	1.500 W
T <sub>max</sub>	1.100 °C	1.200 °C	1.200 °C	1.200 °C	1.200 °C

# DEKEMA. Device equipment.

	AUSTROMAT™ 220	AUSTROMAT™ 624	AUSTROMAT™ 624i™	AUSTROMAT™ 654	AUSTROMAT™ 654i™
Housing: Aluminum					
Corporate Design Customization	–	optional	optional	optional	optional
Touchscreen	4,3", color	7,5", color	10,4", color	7,5", color	10,4", color
USB interface / Ethernet interface	–/–	 / optional		 / optional	
Quartz glass heating element: DEKEMA specific heating element coil guided in pure quartz glass					
Fiber insulation: DEKEMA specific fiber insulation				–	–
DEKEMA stone insulation: internally fabricated, specific stone insulation	–	–	–		
DEKEMA vacuum pump	optional				
DEKEMA ventilation unit	optional	integrated	integrated	integrated	integrated
DEKEMA AutoDry™	–				
MPTA (Multiple Press Time Automatic)	–	–	–	